



Jaundice TL/TL-D

TL-D 18W/52

Tubular or compact single-ended low-pressure mercury vapour fluorescent lamps

Warnings and Safety

- A lamp breaking is extremely unlikely to have any impact on your health. If a lamp breaks, ventilate the room for 30 minutes and remove the parts, preferably with gloves. Put them in a sealed plastic bag and take it to your local waste facilities for recycling. Do not use a vacuum cleaner.

Product data

General Information	
Cap base	G13 [Medium Bi-Pin Fluorescent]
Main application	-
Life to 50% failures (nom.)	13000 h
Useful life (nom.)	2000 h

Light Technical	
Colour Code	52
Lamp Luminous Flux 25°C EL (Nom)	395 lm
Chromaticity coordinate X (nom.)	157
Chromaticity coordinate Y (nom.)	75

Operating and Electrical	
Power (Rated) (Nom)	18 W
Lamp current (nom.)	0.36 A
Voltage (Nom)	59 V

Controls and Dimming	
Dimmable	Yes


Approval and Application	
Mercury (Hg) content (nom.)	13.0 mg

Product Data	
Full product code	872790083485700
Order product name	TL-D 18W/52
EAN/UPC – product	8727900834857
Order code	928048005211
SAP numerator – quantity per pack	1
Numerator – packs per outer box	25
SAP material	928048005211
SAP net weight (piece)	71.000 g

Jaundice TL/TL-D

Dimensional drawing

Product	D (max)	A (max)	B (max)	B (min)	C (max)
TL-D 18W/52	28 mm	589.8 mm	596.9 mm	594.5 mm	604 mm



TL-D 18W/52

